

Logiq C3

Premium







Caractéristiques

Dimensions	50 x 130 x 63 cm
Poids	62 Kg
Ecran	15" LCD

Modes

Configuration standard	Options
Bidimensionnel (B mode)	Logiq View
Harmonique (CHI)	Doppler continu CW (3S RC)
TM (M-Mode)	TM anatomique (AMM)
Doppler pulsé (PW)	

Sondes

Types	Noms	Photos	Fréq. MHz	Applications
µConvexe	6,5C RC		4-10	Echographie abdominale Echocardiographie 2D/TM (0,5-80 Kg)
Phased Array Sonde Doppler	3S-RC		1,5-3,6	Echocardiographie grands animaux
Convexe	C2-5 RC		2-5,5	Echographie abdominale grands animaux
Linéaire	L6-12 RC		4-12	Œil, abdomen petits animaux, NAC, tendons, articulations, thyroïde, vaisseaux périphériques, animaux sauvages, poissons et oiseaux.

Pour plus d'informations sur nos produits, consultez notre site internet

www.scilvet.fr / Tél : 03 90 20 16 40